



808 Laser Hair Removal Machine 12mm*12mm Spotsize With 15 Inch Touch Screen

Our Product Introduction

For more products please visit us on hifubeauty-machine.com


Basic Information

- Place of Origin: China Guangzhou
- Brand Name: 808 hair removal instrument
- Certification: CE
- Model Number: NGY-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: 808 Laser Hair Removal Machine
- Safety Features: Skin Contact Cooling, Automatic Power Cut-off
- Spotsize: 12mm*12mm
- Skin Types: I-VI
- Machine App: Permanent Laser Hair Removal Machines
- Display: 15-inch Touch Screen
- Power: 1000W
- Weight: 15 Kg
- Pulse Duration: 10-400 Ms
- Highlight: 808 Laser Hair Removal Machine,
808nm Laser Hair Removal Machine,
Laser Hair Removal Machine 808nm



More Images





Product Description

Advanced 808 Laser Hair Removal Machine With 15-Inch Touch Screen

Product Description:

The device boasts a user-friendly Display: 15-inch Touch Screen, which makes it easy to operate. You can customize the settings according to your hair type and skin tone. The device's Pulse Duration: 10-400 Ms feature ensures that the laser targets the hair follicles effectively, providing long-lasting results.

With the Package Size: 75*54*54cm³, this machine is compact and easy to store. You can use it in the comfort of your home or take it with you on your travels.

The machine's Power: 1200W is powerful enough to remove even the thickest hair. The device's Safety Features: Skin Contact Cooling, Automatic Power Cut-off, ensures that the machine is safe to use. The skin contact cooling feature protects your skin from burns, while the automatic power cut-off feature prevents the device from overheating.







Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou